



C-Bona Closed Suction System C-Bona 72 Closed Suction System

The closed suction system is used with patients who are mechanically ventilated and have an endotracheal or tracheostomy tube in place. The advantage of the closed method of suctioning is that it doesn't require disconnecting the ventilator or oxygen source from the artificial airway when patient needs to be suctioned. This reduces the risks of infection and oxygen desaturation.

C-Bona Closed Suction System

- Equipped with irrigation port with one way valve
- Equipped with MDI port
- Equipped with double swivel connector
- Available for 24hrs and 72hrs of continuous use.



ISO color coded for instant recognition of catheter sizes

Once the airway is clear, withdraw the catheter completely into the plastic sheath. A black mark on the catheter indicates full withdrawal.



With the unique isolation valve, C-Bona 72 Closed Suction System features safe and simple access for mechanical ventilated patients for 72 hours of continuous use.

Smooth tubing surface and atraumatic rounded catheter tip for minimal tracheal wall injury.

C-Bona Closed Suction System

Product Code	Outer Diameter	Cath. length	Remark	Packaging	ISO Color
CB-6001000A / CB-6001000AM CB-6021000A / CB-6021000AM CB-6011000A / CB-6011000AM CB-6031000A / CB-6031000AM	10FR / 3.33mm	540mm	Endotracheal	20ea/box	●
CB-6001200 / CB-6001200M CB-6021200 / CB-6021200M CB-6011200 / CB-6011200M CB-6031200 / CB-6031200M		305mm	Tracheostomy		
CB-6001400 / CB-6001400M CB-6021400 / CB-6021400M CB-6011400 / CB-6011400M CB-6031400 / CB-6031400M	12FR / 4.0mm	540mm	Endotracheal	20ea/box	○
CB-6001600 / CB-6001600M CB-6021600 / CB-6021600M CB-6011600 / CB-6011600M CB-6031600 / CB-6031600M		305mm	Tracheostomy		
CB-6001800 / CB-6001800M CB-6021800 / CB-6021800M CB-6011800 / CB-6011800M CB-6031800 / CB-6031800M	14FR / 4.67mm	540mm	Endotracheal	20ea/box	●
CB-6002000 / CB-6002000M CB-6022000 / CB-6022000M CB-6012000 / CB-6012000M CB-6032000 / CB-6032000M		305mm	Tracheostomy		
CB-6002200 / CB-6002200M CB-6022200 / CB-6022200M CB-6012200 / CB-6012200M CB-6032200 / CB-6032200M	16FR / 5.3mm	540mm	Endotracheal	20ea/box	●
CB-6002400 / CB-6002400M CB-6022400 / CB-6022400M CB-6012400 / CB-6012400M CB-6032400 / CB-6032400M		305mm	Tracheostomy		

C-Bona 72 Closed Suction System

Product Code	Outer Diameter	Cath. length	Remark	Packaging	ISO Color
CB-6051000A / CB-6051000AM CB-6071000A / CB-6071000AM CB-6061000A / CB-6061000AM CB-6081000A / CB-6081000AM	10FR / 3.33mm	540mm	Endotracheal	20ea/box	●
CB-6051200 / CB-6051200M CB-6071200 / CB-6071200M CB-6061200 / CB-6061200M CB-6081200 / CB-6081200M		305mm	Tracheostomy		
CB-6051400 / CB-6051400M CB-6071400 / CB-6071400M CB-6061400 / CB-6061400M CB-6081400 / CB-6081400M	12FR / 4.0mm	540mm	Endotracheal	20ea/box	○
CB-6051600 / CB-6051600M CB-6071600 / CB-6071600M CB-6061600 / CB-6061600M CB-6081600 / CB-6081600M		305mm	Tracheostomy		
CB-6051800 / CB-6051800M CB-6071800 / CB-6071800M CB-6061800 / CB-6061800M CB-6081800 / CB-6081800M	14FR / 4.67mm	540mm	Endotracheal	20ea/box	●
CB-6052000 / CB-6052000M CB-6072000 / CB-6072000M CB-6062000 / CB-6062000M CB-6082000 / CB-6082000M		305mm	Tracheostomy		
CB-6052200 / CB-6052200M CB-6072200 / CB-6072200M CB-6062200 / CB-6062200M CB-6082200 / CB-6082200M	16FR / 5.3mm	540mm	Endotracheal	20ea/box	●
CB-6052400 / CB-6052400M CB-6072400 / CB-6072400M CB-6062400 / CB-6062400M CB-6082400 / CB-6082400M		305mm	Tracheostomy		

■ M = with MDI port



Vitaltec Vitaltec Corporation

No. 12, Lane 4-30, Chyuan-Zhou Rd., Hou-Li Dist., 42142 Taichung City, TAIWAN

Website: www.vitaltec.com.tw
TEL: +886-4-2558-0886

Made in Taiwan

